<u>Liquid Media / Swabs to be used for South of Tyne Microbiology Tests</u>



Pink Capped Tops
(with additional swabs for multiple sites) mainly used for MRSA screening swabs Red swabs must be discarded after application in liquid.

application in liquid. 2 or 3 swab variants available depending upon the requirement



Orange Capped
Tops for General
Microbiology
Bacterial/Fungal
Culture. Mini tip
swab for Urethral or
ENT use only



Pink Capped Tops
(most commonly
used) for General
Microbiology
Bacterial/Fungal
Culture. Single
swab for individual
sites



Green Capped
Tops for Viral
Investigations ONLY
(incl COVID19 and
Flu). These cannot
be used for
bacterial culture



Yellow Capped
Tops - cobas media
for Chlamydia
Gonorrhoea,
Trichomonas or
Mycoplasma PCR
ONLY





Existing Blue
topped Gel swabs
are to be phased
out across the
South of Tyne and
Wear Hospital
and Primary Care
sites from
February 2025

